

T-cell surface antigen CD2

Catalog No: #AP81599



Package Size: #AP81599-1 50ug #AP81599-2 100ug #AP81599-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface antigen CD2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LFA-2,LFA-3 receptor,Lymphocyte antigen 37,Ly-37,T-cell surface antigen T11/Leu-5,CD2,Ly-37
Accession No.	Uniprot ID:P08920
Uniprot	P08920
GeneID	12481;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd2

Note: This product is for in vitro research use only