

Calcitonin gene-related peptide type 1 receptor

Catalog No: #AP81601



Package Size: #AP81601-1 50ug #AP81601-2 100ug #AP81601-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcitonin gene-related peptide type 1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CGRP type 1 receptor, Calcitonin receptor-like receptor, CGRPR
Accession No.	Uniprot ID:Q8WN93
Uniprot	Q8WN93
GeneID	397172;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CALCRL

Note: This product is for in vitro research use only