

## Caprin-1

Catalog No: #AP81604



Package Size: #AP81604-1 50ug #AP81604-2 100ug #AP81604-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Caprin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell cycle-associated protein 1,Cytoplasmic activation- and proliferation-associated protein 1,GPI-anchored membrane protein 1,GPI-anchored protein p137,GPI-p137,Membrane component chromosome 11 surface marker 1,RNA granule protein 105,GPIAP1,GPIP137,M11S1,RNG105
Accession No.	Uniprot ID:Q14444
Uniprot	Q14444
GeneID	4076;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CAPRIN1

Note: This product is for in vitro research use only