

Dipeptidyl peptidase 1

Catalog No: #AP81605



Package Size: #AP81605-1 50ug #AP81605-2 100ug #AP81605-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Dipeptidyl peptidase 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Cathepsin C,Cathepsin J,Dipeptidyl peptidase I,DPP-I,Dipeptidyl transferase |
| Accession No. | Uniprot ID:Q3ZCJ8 |
| Uniprot | Q3ZCJ8 |
| GeneID | 352958; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTSC

Note: This product is for in vitro research use only