

# Cell division cycle protein 23 homolog

Catalog No: #AP81606



Package Size: #AP81606-1 50ug #AP81606-2 100ug #AP81606-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cell division cycle protein 23 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Anaphase-promoting complex subunit 8,APC8,Cyclosome subunit 8,Anapc8
Accession No.	Uniprot ID:Q8BGZ4
Uniprot	Q8BGZ4
GeneID	52563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cdc23

Note: This product is for in vitro research use only