

Homeobox protein CDX-2

Catalog No: #AP81607



Package Size: #AP81607-1 50ug #AP81607-2 100ug #AP81607-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Homeobox protein CDX-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Caudal-type homeobox protein 2, Cdx-2
Accession No.	Uniprot ID:P43241
Uniprot	P43241
GeneID	12591;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cdx2

Note: This product is for in vitro research use only