

# Eotaxin

Catalog No: #AP81611



Package Size: #AP81611-1 50ug #AP81611-2 100ug #AP81611-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Eotaxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-C motif chemokine 11,Eosinophil chemotactic protein,Small-inducible cytokine A11,Scya11
Accession No.	Uniprot ID:P48298
Uniprot	P48298
GeneID	20292;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ccl11

Note: This product is for in vitro research use only