

CD2-associated protein

Catalog No: #AP81613



Package Size: #AP81613-1 50ug #AP81613-2 100ug #AP81613-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD2-associated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Adapter protein CMS,Cas ligand with multiple SH3 domains
Accession No.	Uniprot ID:Q9Y5K6
Uniprot	Q9Y5K6
GeneID	23607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD2AP

Note: This product is for in vitro research use only