

60 kDa heat shock protein, mitochondrial

Catalog No: #AP81616



Package Size: #AP81616-1 50ug #AP81616-2 100ug #AP81616-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	60 kDa heat shock protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	60 kDa chaperonin,Chaperonin 60,CPN60,Heat shock protein 60,HSP-60,HuCHA60,Mitochondrial matrix protein P1,P60 lymphocyte protein,HSP60
Accession No.	Uniprot ID:P10809
Uniprot	P10809
GeneID	3329;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPD1

Note: This product is for in vitro research use only