

# CD70 antigen

Catalog No: #AP81619



Package Size: #AP81619-1 50ug #AP81619-2 100ug #AP81619-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD70 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CD27 ligand,CD27-L,Tumor necrosis factor ligand superfamily member 7,CD70,CD27L,TNFSF7
Accession No.	Uniprot ID:Q3ZDR4
Uniprot	Q3ZDR4
GeneID	733577;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD70

Note: This product is for in vitro research use only