

## Signal transducer CD24

Catalog No: #AP81620



Package Size: #AP81620-1 50ug #AP81620-2 100ug #AP81620-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Signal transducer CD24
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Small cell lung carcinoma cluster 4 antigen,CD24,CD24A
Accession No.	Uniprot ID:P25063
Uniprot	P25063
GeneID	100133941;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD24

Note: This product is for in vitro research use only