

Chitinase-3-like protein 1

Catalog No: #AP81622



Package Size: #AP81622-1 50ug #AP81622-2 100ug #AP81622-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chitinase-3-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BRP39 protein,Breast regression protein 39,Cartilage glycoprotein 39,CGP-39,Brp39,Chil1
Accession No.	Uniprot ID:Q61362
Uniprot	Q61362
GeneID	12654;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Chi311

Note: This product is for in vitro research use only