

C-C motif chemokine 25

Catalog No: #AP81625



Package Size: #AP81625-1 50ug #AP81625-2 100ug #AP81625-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------------------------------------------------------------|
| Product Name | C-C motif chemokine 25 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Chemokine TECK, Small-inducible cytokine A25, Thymus-expressed chemokine, Scya25, Teck |
| Accession No. | Uniprot ID:O35903 |
| Uniprot | O35903 |
| GeneID | 20300; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ccl25

Note: This product is for in vitro research use only