

## CD40 ligand

Catalog No: #AP81626



Package Size: #AP81626-1 50ug #AP81626-2 100ug #AP81626-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD40 ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD40-L,T-cell antigen Gp39,TNF-related activation protein,TRAP,Tumor necrosis factor ligand superfamily member 5,CD154,CD40L,TNFSF5,TRAP
Accession No.	Uniprot ID:P29965
Uniprot	P29965
GeneID	959;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD40LG

Note: This product is for in vitro research use only