

Chymotrypsin-like elastase family member 1

Catalog No: #AP81635

Package Size: #AP81635-1 50ug #AP81635-2 100ug #AP81635-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chymotrypsin-like elastase family member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Elastase I,Elastase-1,ELA1
Accession No.	Uniprot ID:Q28153
Uniprot	Q28153
GeneID	281139;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CELA1

Note: This product is for in vitro research use only