

C-C motif chemokine 22

Catalog No: #AP81638



Package Size: #AP81638-1 50ug #AP81638-2 100ug #AP81638-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 22
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CC chemokine STCP-1,MDC(1-69),Macrophage-derived chemokine,Small-inducible cytokine A22,Stimulated T-cell chemotactic protein 1,MDC,SCYA22,A-152E5.1
Accession No.	Uniprot ID:O00626
Uniprot	O00626
GeneID	6367;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL22

Note: This product is for in vitro research use only