

# CD209 antigen

Catalog No: #AP81640



Package Size: #AP81640-1 50ug #AP81640-2 100ug #AP81640-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD209 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-type lectin domain family 4 member L,Dendritic cell-specific ICAM-3-grabbing non-integrin 1,DC-SIGN,CD209,CLEC4L
Accession No.	Uniprot ID:Q9NNX6
Uniprot	Q9NNX6
GeneID	30835;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD209

Note: This product is for in vitro research use only