

# T-cell antigen CD7

Catalog No: #AP81641



Package Size: #AP81641-1 50ug #AP81641-2 100ug #AP81641-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	T-cell antigen CD7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GP40,T-cell leukemia antigen,T-cell surface antigen Leu-9,TP41,CD7
Accession No.	Uniprot ID:P09564
Uniprot	P09564
GeneID	924;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD7

Note: This product is for in vitro research use only