

Cyclin-dependent kinase inhibitor 1

Catalog No: #AP81644



Package Size: #AP81644-1 50ug #AP81644-2 100ug #AP81644-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase inhibitor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CDK-interacting protein 1,Melanoma differentiation-associated protein,p21,Cip1,Waf1
Accession No.	Uniprot ID:P39689
Uniprot	P39689
GeneID	12575;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdkn1a

Note: This product is for in vitro research use only