

C-type lectin domain family 4 member E

Catalog No: #AP81646



Package Size: #AP81646-1 50ug #AP81646-2 100ug #AP81646-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | C-type lectin domain family 4 member E |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | C-type lectin superfamily member 9, Macrophage-inducible C-type lectin, CLECSF9, MINCLE, UNQ218/PRO244 |
| Accession No. | Uniprot ID:Q9ULY5 |
| Uniprot | Q9ULY5 |
| GeneID | 26253; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CLEC4E

Note: This product is for in vitro research use only