

Collectin-11

Catalog No: #AP81648



Package Size: #AP81648-1 50ug #AP81648-2 100ug #AP81648-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Collectin-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Collectin kidney protein 1,CL-K1
Accession No.	Uniprot ID:Q3SXB8
Uniprot	Q3SXB8
GeneID	71693;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Colec11

Note: This product is for in vitro research use only