

TTR Mouse Monoclonal Antibody(1D7)

Catalog No: #44227

Package Size: #44227-1 50ul #44227-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	TTR Mouse Monoclonal Antibody(1D7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1D7
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous TTR protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein of TTR
Target Name	TTR
Other Names	Transthyretin (ATTR) (Prealbumin) (TBPA)
Accession No.	Swiss-Prot:P02766Gene ID:7276
Uniprot	P02766
GeneID	7276;
Calculated MW	16kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:500-1:2000

Immunohistochemistry: 1:50-1:2000


Images



Western blot analysis of Human Serum using TTR Mouse mAb diluted at 1:2000



Immunohistochemical analysis of paraffin-embedded Human Liver Carcinoma Tissue using TTR Mouse mAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Thyroid Tissue using TTR Mouse mAb diluted at 1:200

Note: This product is for in vitro research use only