

C-C chemokine receptor type 1

Catalog No: #AP81654



Package Size: #AP81654-1 50ug #AP81654-2 100ug #AP81654-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-C CKR-1,Macrophage inflammatory protein 1-alpha receptor,MIP-1alpha-R,RANTES-R,CD191,Cmkbr1
Accession No.	Uniprot ID:P51675
Uniprot	P51675
GeneID	12768;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccr1

Note: This product is for in vitro research use only