

C-C chemokine receptor type 5

Catalog No: #AP81655



Package Size: #AP81655-1 50ug #AP81655-2 100ug #AP81655-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	C-C CKR-5,CC-CKR-5,CCR-5,CCR5,CMKBR5,CD195
Accession No.	Uniprot ID:Q5ECR9
Uniprot	Q5ECR9
GeneID	484789;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCR5

Note: This product is for in vitro research use only