

C-C motif chemokine 7

Catalog No: #AP81657



Package Size: #AP81657-1 50ug #AP81657-2 100ug #AP81657-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Monocyte chemoattractant protein 3, Monocyte chemotactic protein 3, MCP-3, NC28, Small-inducible cytokine A7, MCP3, SCYA6, SCYA7
Accession No.	Uniprot ID:P80098
Uniprot	P80098
GeneID	6354;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL7

Note: This product is for in vitro research use only