

Monocyte differentiation antigen CD14

Catalog No: #AP81659

Package Size: #AP81659-1 50ug #AP81659-2 100ug #AP81659-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Monocyte differentiation antigen CD14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Myeloid cell-specific leucine-rich glycoprotein,CD14
Accession No.	Uniprot ID:Q95122
Uniprot	Q95122
GeneID	281048;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD14

Note: This product is for in vitro research use only