

TTR Mouse Monoclonal Antibody(2B12)

Catalog No: #44228

Package Size: #44228-1 50ul #44228-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|-----------------------|---|
| Product Name | TTR Mouse Monoclonal Antibody(2B12) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 2B12 |
| Isotype | IgG1 |
| Applications | WB IHC |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous TTR protein |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant Protein of TTR |
| Target Name | TTR |
| Other Names | Transthyretin (ATTR) (Prealbumin) (TBPA) |
| Accession No. | Swiss-Prot:P02766Gene ID:7276 |
| Uniprot | P02766 |
| GeneID | 7276; |
| Calculated MW | 16kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details


Western blotting: 1:500-1:2000

Immunohistochemistry: 1:50-1:2000

Images



Western blot analysis of Human Serum using TTR Mouse mAb diluted at 1:2000



Immunohistochemical analysis of paraffin-embedded Human Thyroid Tissue using TTR Mouse mAb diluted at 1:200

Note: This product is for in vitro research use only