

Monocyte differentiation antigen CD14

Catalog No: #AP81660



Package Size: #AP81660-1 50ug #AP81660-2 100ug #AP81660-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Monocyte differentiation antigen CD14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Myeloid cell-specific leucine-rich glycoprotein,CD14
Accession No.	Uniprot ID:Q63691
Uniprot	Q63691
GeneID	60350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd14

Note: This product is for in vitro research use only