

CD166 antigen

Catalog No: #AP81661



Package Size: #AP81661-1 50ug #AP81661-2 100ug #AP81661-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | CD166 antigen |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Activated leukocyte cell adhesion molecule, SB-10 antigen, CD166 |
| Accession No. | Uniprot ID:O46634 |
| Uniprot | O46634 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ALCAM

Note: This product is for in vitro research use only