

## Cytochrome P450 1A1

Catalog No: #AP81664



Package Size: #AP81664-1 50ug #AP81664-2 100ug #AP81664-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytochrome P450 1A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CYP1A1,Cytochrome P450 form 6,Cytochrome P450-C,Cytochrome P450-P1
Accession No.	Uniprot ID:P04798
Uniprot	P04798
GeneID	1543;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CYP1A1

Note: This product is for in vitro research use only