

Caspase-2

Catalog No: #AP81667



Package Size: #AP81667-1 50ug #AP81667-2 100ug #AP81667-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CASP-2,Neural precursor cell expressed developmentally down-regulated protein 2,NEDD-2,Protease ICH-1,Ich1,Nedd-2,Nedd2
Accession No.	Uniprot ID:P29594
Uniprot	P29594
GeneID	12366;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Casp2

Note: This product is for in vitro research use only