

Bile salt-activated lipase

Catalog No: #AP81670



Package Size: #AP81670-1 50ug #AP81670-2 100ug #AP81670-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bile salt-activated lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	BAL, Bile salt-stimulated lipase, BSSL, Carboxyl ester lipase, Cholesterol esterase, Pancreatic lysophospholipase, Sterol esterase
Accession No.	Uniprot ID: P07882
Uniprot	P07882
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cel

Note: This product is for in vitro research use only