

## CCAAT/enhancer-binding protein beta

Catalog No: #AP81671



Package Size: #AP81671-1 50ug #AP81671-2 100ug #AP81671-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CCAAT/enhancer-binding protein beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C/EBP beta,Liver activator protein,LAP,Liver-enriched inhibitory protein,LIP,Nuclear factor NF-IL6,Transcription factor 5,TCF-5,TCF5,PP9092
Accession No.	Uniprot ID:P17676
Uniprot	P17676
GeneID	1051;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CEBPB

Note: This product is for in vitro research use only