

Signal transducer CD24

Catalog No: #AP81673



Package Size: #AP81673-1 50ug #AP81673-2 100ug #AP81673-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Signal transducer CD24
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lymphocyte antigen 52, Ly-52, M1/69-J11D heat stable antigen, HSA, Nectadrin, R13-Ag, X62 heat stable antigen, CD24, Cd24a, Ly-52
Accession No.	Uniprot ID: P24807
Uniprot	P24807
GeneID	12484;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd24

Note: This product is for in vitro research use only