

Endosialin

Catalog No: #AP81674



Package Size: #AP81674-1 50ug #AP81674-2 100ug #AP81674-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endosialin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Tumor endothelial marker 1,CD248,CD164L1,TEM1
Accession No.	Uniprot ID:Q9HCU0
Uniprot	Q9HCU0
GeneID	57124;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD248

Note: This product is for in vitro research use only