

CD70 antigen

Catalog No: #AP81675



Package Size: #AP81675-1 50ug #AP81675-2 100ug #AP81675-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD70 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD27 ligand,CD27-L,Tumor necrosis factor ligand superfamily member 7,CD70,CD27L,CD27LG,TNFSF7
Accession No.	Uniprot ID:P32970
Uniprot	P32970
GeneID	970;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD70

Note: This product is for in vitro research use only