

Cyclin-dependent kinase 2

Catalog No: #AP81676



Package Size: #AP81676-1 50ug #AP81676-2 100ug #AP81676-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cell division protein kinase 2,CDKN2
Accession No.	Uniprot ID:Q5E9Y0
Uniprot	Q5E9Y0
GeneID	519217;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CDK2

Note: This product is for in vitro research use only