

Carbonic anhydrase 4

Catalog No: #AP81679



Package Size: #AP81679-1 50ug #AP81679-2 100ug #AP81679-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carbonic anhydrase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carbonate dehydratase IV,Carbonic anhydrase IV,CA-IV,Car4
Accession No.	Uniprot ID:Q64444
Uniprot	Q64444
GeneID	12351;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ca4

Note: This product is for in vitro research use only