

Ceramide kinase

Catalog No: #AP81681



Package Size: #AP81681-1 50ug #AP81681-2 100ug #AP81681-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ceramide kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mCERK,Acylsphingosine kinase
Accession No.	Uniprot ID:Q8K4Q7
Uniprot	Q8K4Q7
GeneID	223753;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cerk

Note: This product is for in vitro research use only