

T-cell surface glycoprotein CD3 gamma chain

Catalog No: #AP81682



Package Size: #AP81682-1 50ug #AP81682-2 100ug #AP81682-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD3 gamma chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	T-cell receptor T3 gamma chain,CD3g
Accession No.	Uniprot ID:Q5PXD3
Uniprot	Q5PXD3
GeneID	494013;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD3G

Note: This product is for in vitro research use only