

CD40 ligand

Catalog No: #AP81683



Package Size: #AP81683-1 50ug #AP81683-2 100ug #AP81683-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD40 ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD154, CD40L, TNFSF5
Accession No.	Uniprot ID:Q95MQ5
Uniprot	Q95MQ5
GeneID	397231;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CD40LG

Note: This product is for in vitro research use only