

Cerebral dopamine neurotrophic factor

Catalog No: #AP81686



Package Size: #AP81686-1 50ug #AP81686-2 100ug #AP81686-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cerebral dopamine neurotrophic factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ARMET-like protein 1,Arginine-rich protein mutated in early stage tumors-like 1,Conserved dopamine neurotrophic factor,Armetl1
Accession No.	Uniprot ID:P0C5I0
Uniprot	P0C5I0
GeneID	361276;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdnf

Note: This product is for in vitro research use only