

Cytidine deaminase

Catalog No: #AP81687



Package Size: #AP81687-1 50ug #AP81687-2 100ug #AP81687-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytidine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytidine aminohydrolase,CDD
Accession No.	Uniprot ID:P32320
Uniprot	P32320
GeneID	978;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDA

Note: This product is for in vitro research use only