

## Cyclin-dependent kinase 12

Catalog No: #AP81691



Package Size: #AP81691-1 50ug #AP81691-2 100ug #AP81691-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cyclin-dependent kinase 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cdc2-related kinase, arginine/serine-rich,CrkRS,Cell division cycle 2-related protein kinase 7,CDC2-related protein kinase 7,Cell division protein kinase 12,Protein kinase for splicing component,Crk7,Crkrs,Pksc
Accession No.	Uniprot ID:Q3MJK5
Uniprot	Q3MJK5
GeneID	192350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cdk12

Note: This product is for in vitro research use only