

CAMPATH-1 antigen

Catalog No: #AP81694



Package Size: #AP81694-1 50ug #AP81694-2 100ug #AP81694-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CAMPATH-1 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lymphocyte differentiation antigen B7,CD52,Cdw52,Mb7
Accession No.	Uniprot ID:Q64389
Uniprot	Q64389
GeneID	23833;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd52

Note: This product is for in vitro research use only