

CD59A glycoprotein

Catalog No: #AP81695



Package Size: #AP81695-1 50ug #AP81695-2 100ug #AP81695-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD59A glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAC-inhibitory protein,MAC-IP,Membrane attack complex inhibition factor,MACIF,Protectin,CD59,Cd59
Accession No.	Uniprot ID:O55186
Uniprot	O55186
GeneID	12509;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd59a

Note: This product is for in vitro research use only