

Bile salt-activated lipase

Catalog No: #AP81702



Package Size: #AP81702-1 50ug #AP81702-2 100ug #AP81702-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bile salt-activated lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	BAL,Bile salt-stimulated lipase,BSSL,Carboxyl ester lipase,Cholesterol esterase,Pancreatic lysophospholipase,Sterol esterase
Accession No.	Uniprot ID:P30122
Uniprot	P30122
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CEL

Note: This product is for in vitro research use only