

Eotaxin

Catalog No: #AP81704



Package Size: #AP81704-1 50ug #AP81704-2 100ug #AP81704-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Eotaxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C-C motif chemokine 11,Eosinophil chemotactic protein,Small-inducible cytokine A11,Scya11
Accession No.	Uniprot ID:P97545
Uniprot	P97545
GeneID	29397;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccl11

Note: This product is for in vitro research use only