

C-C motif chemokine 8

Catalog No: #AP81706

Package Size: #AP81706-1 50ug #AP81706-2 100ug #AP81706-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Monocyte chemoattractant protein 2, Monocyte chemotactic protein 2, MCP-2, Small-inducible cytokine A8, Mcp2, Scya8
Accession No.	Uniprot ID:Q9Z121
Uniprot	Q9Z121
GeneID	20307;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ccl8

Note: This product is for in vitro research use only