

Cyclin-dependent kinase 16

Catalog No: #AP81709



Package Size: #AP81709-1 50ug #AP81709-2 100ug #AP81709-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell division protein kinase 16,PCTAIRE-motif protein kinase 1,Serine/threonine-protein kinase PCTAIRE-1,PCTAIRE1,PCTK1
Accession No.	Uniprot ID:Q00536
Uniprot	Q00536
GeneID	5127;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDK16

Note: This product is for in vitro research use only